

(R)-1,2-Distearin

Chemical Properties


CAS No. : 1429-59-0

Formula: C₃₉H₇₆O₅

Molecular Weight: 625.02

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R)-1,2-Distearin is a metabolic intermediate of natural triglycerides and a stereoisomer of synthetically produced diacylglycerol (DAG). (R)-1,2-Distearin is a nonionic, fat-soluble small molecule.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5999 mL	7.9997 mL	15.9995 mL
5 mM	0.320 mL	1.5999 mL	3.1999 mL
10 mM	0.160 mL	0.800 mL	1.5999 mL
50 mM	0.032 mL	0.160 mL	0.320 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481