

Besilesomab

Chemical Properties

CAS No. : 537694-98-7

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Besilesomab is a monoclonal antibody targeting CEACAM8 (CD66b) and related carcinoembryonic antigen family members. Besilesomab can be used for neutrophil biology and oncology research.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481