

## Anti-Mouse IFNAR-1 Antibody (15A3.B)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse IFNAR-1 Antibody (15A3.B) is a monoclonal antibody targeting mouse IFNAR-1. It can be used for immunology research involving interferon signaling.
-------------	---

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481