

## Anti-Mouse IL-17A Antibody (3A7.N)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Store at 4°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Mouse IL-17A Antibody (3A7.N) is a monoclonal antibody targeting mouse IL-17A. It can be used for Th17 pathway and inflammation research.
-------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481