

## Human IgG1 kappa-MMAE (PEG8), ADC Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Store at 4°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human IgG1 kappa-MMAE (PEG8), ADC Control is an IgG1-based antibody-drug conjugate control containing MMAE (PEG8). It can be used as a control in ADC.
-------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481