

## Dodec-2E-ene-8,10-diynoic acid isobutylamide

## Chemical Properties


CAS No. : 120727-29-9

Formula: C<sub>16</sub>H<sub>23</sub>NO

Molecular Weight: 245.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dodec-2E-ene-8,10-diynoic acid isobutylamide is a natural product that can be used in life science research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0756 mL	20.3782 mL	40.7564 mL
5 mM	0.8151 mL	4.0756 mL	8.1513 mL
10 mM	0.4076 mL	2.0378 mL	4.0756 mL
50 mM	0.0815 mL	0.4076 mL	0.8151 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481