

Methyl glycyrrhetinate

Chemical Properties


CAS No. : 1477-44-7

Formula: C₃₁H₄₈O₄

Molecular Weight: 484.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl glycyrrhetinate is a biologically active compound that can be utilized for life science research applications.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0631 mL	10.3154 mL	20.6309 mL
5 mM	0.4126 mL	2.0631 mL	4.1262 mL
10 mM	0.2063 mL	1.0315 mL	2.0631 mL
50 mM	0.0413 mL	0.2063 mL	0.4126 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481