


3'-p-Hydroxypaclitaxel

Chemical Properties

CAS No. :	132160-32-8
Formula:	C47H51NO15
Molecular Weight:	869.91
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-p-Hydroxypaclitaxel is a natural product that can be used in life science research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1495 mL	5.7477 mL	11.4954 mL
5 mM	0.2299 mL	1.1495 mL	2.2991 mL
10 mM	0.115 mL	0.5748 mL	1.1495 mL
50 mM	0.023 mL	0.115 mL	0.2299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481