


## Isorebaudioside A

## Chemical Properties

CAS No. :	1058600-85-3
Formula:	C <sub>44</sub> H <sub>70</sub> O <sub>23</sub>
Molecular Weight:	967.01
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Isorebaudioside A is a commonly used biochemical reagent that can serve as a biomaterial or organic compound for life science research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0341 mL	5.1706 mL	10.3412 mL
5 mM	0.2068 mL	1.0341 mL	2.0682 mL
10 mM	0.1034 mL	0.5171 mL	1.0341 mL
50 mM	0.0207 mL	0.1034 mL	0.2068 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481