

Compound T129299

Chemical Properties

CAS No. :

Formula: C₃₀H₅₄O₁₆

Molecular Weight: 670.75

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	16-(((2R,3R,4R,5S,6R)-6-[(acetyloxy)methyl]-3,4-dihydroxy-5-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)oxan-2-yl]oxy)-3,15-dihydroxyhexadecanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4909 mL	7.4543 mL	14.9087 mL
5 mM	0.2982 mL	1.4909 mL	2.9817 mL
10 mM	0.1491 mL	0.7454 mL	1.4909 mL
50 mM	0.0298 mL	0.1491 mL	0.2982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481