


Fumarprotocetraric acid

Chemical Properties

CAS No. :	489-50-9
Formula:	C ₂₂ H ₁₆ O ₁₂
Molecular Weight:	472.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fumarprotocetraric acid is a natural product that can be used in life science research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.117 mL	10.5851 mL	21.1703 mL
5 mM	0.4234 mL	2.117 mL	4.2341 mL
10 mM	0.2117 mL	1.0585 mL	2.117 mL
50 mM	0.0423 mL	0.2117 mL	0.4234 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481