


Sempervirine nitric acid

Chemical Properties

CAS No. :	5436-46-4
Formula:	C ₁₉ H ₁₇ N ₃ O ₃
Molecular Weight:	335.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sempervirine Nitrate is a commonly used biochemical reagent that can serve as a biomaterial or organic compound for life science research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9819 mL	14.9094 mL	29.8187 mL
5 mM	0.5964 mL	2.9819 mL	5.9637 mL
10 mM	0.2982 mL	1.4909 mL	2.9819 mL
50 mM	0.0596 mL	0.2982 mL	0.5964 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481