


2H-1,3,5-thiadiazine-2-thione,tetrahydro-3,5-dipropyl-

Chemical Properties

CAS No. :	52416-62-3
Formula:	C ₉ H ₁₈ N ₂ S ₂
Molecular Weight:	218.39
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2H-1,3,5-thiadiazine-2-thione,tetrahydro-3,5-dipropyl- ,with CAS No. 52416-62-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2H-1,3,5-thiadiazine-2-thione,tetrahydro-3,5-dipropyl- provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.579 mL	22.8948 mL	45.7896 mL
5 mM	0.9158 mL	4.579 mL	9.1579 mL
10 mM	0.4579 mL	2.2895 mL	4.579 mL
50 mM	0.0916 mL	0.4579 mL	0.9158 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481