

2-(4,7-dimethyl-1,2,4a,5,6,8a-hexahydronaphthalen-1-yl)prop-2-enoic acid

Chemical Properties


CAS No. :

Formula: C₁₅H₂₀O₂

Molecular Weight: 232.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(4,7-dimethyl-1,2,4a,5,6,8a-hexahydronaphthalen-1-yl)prop-2-enoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3044 mL	21.522 mL	43.0441 mL
5 mM	0.8609 mL	4.3044 mL	8.6088 mL
10 mM	0.4304 mL	2.1522 mL	4.3044 mL
50 mM	0.0861 mL	0.4304 mL	0.8609 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481