

(4aS,6S)-4-(hydroxymethyl)-3,4a,8,8-tetramethyl-1,4,4a,5,6,7,8,8a-octahydronaphthalene-1,4,6-triol


## Chemical Properties

CAS No. :

Formula: C<sub>15</sub>H<sub>26</sub>O<sub>4</sub>

Molecular Weight: 270.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(4aS,6S)-4-(hydroxymethyl)-3,4a,8,8-tetramethyl-1,4,4a,5,6,7,8,8a-octahydronaphthalene-1,4,6-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6986 mL	18.4932 mL	36.9864 mL
5 mM	0.7397 mL	3.6986 mL	7.3973 mL
10 mM	0.3699 mL	1.8493 mL	3.6986 mL
50 mM	0.074 mL	0.3699 mL	0.7397 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481