

(3S,6S)-3-methyl-6-[[7-(3-methylbut-2-en-1-yl)-1H-indol-3-yl]methyl]piperazine-2,5-dione

Chemical Properties


CAS No. :

Formula: C₁₉H₂₃N₃O₂

Molecular Weight: 325.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(3S,6S)-3-methyl-6-[[7-(3-methylbut-2-en-1-yl)-1H-indol-3-yl]methyl]piperazine-2,5-dione is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.073 mL	15.3652 mL	30.7305 mL
5 mM	0.6146 mL	3.073 mL	6.1461 mL
10 mM	0.3073 mL	1.5365 mL	3.073 mL
50 mM	0.0615 mL	0.3073 mL	0.6146 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481