

## 8-hydroxy-3-(2-hydroxypropyl)-6-methoxy-2H-chromen-2-one

### Chemical Properties


CAS No. :

Formula: C<sub>13</sub>H<sub>14</sub>O<sub>5</sub>

Molecular Weight: 250.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	8-hydroxy-3-(2-hydroxypropyl)-6-methoxy-2H-chromen-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.996 mL	19.980 mL	39.960 mL
5 mM	0.7992 mL	3.996 mL	7.992 mL
10 mM	0.3996 mL	1.998 mL	3.996 mL
50 mM	0.0799 mL	0.3996 mL	0.7992 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481