

## 3-[(3-[(2R,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)decanoyl]oxy]dodecanoic acid

## Chemical Properties


CAS No. :

Formula: C<sub>28</sub>H<sub>52</sub>O<sub>9</sub>

Molecular Weight: 532.72

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-[(3-[(2R,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)decanoyl]oxy]dodecanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8772 mL	9.3858 mL	18.7716 mL
5 mM	0.3754 mL	1.8772 mL	3.7543 mL
10 mM	0.1877 mL	0.9386 mL	1.8772 mL
50 mM	0.0375 mL	0.1877 mL	0.3754 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481