

3-({3-[(3-hydroxydecanoyl)oxy]dodecanoyl}oxy)decanoic acid

Chemical Properties


CAS No. :

Formula: C₃₂H₆₀O₇

Molecular Weight: 556.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|-------------|---|
| Description | 3-({3-[(3-hydroxydecanoyl)oxy]dodecanoyl}oxy)decanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research. |
|-------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.7959 mL | 8.9794 mL | 17.9588 mL |
| 5 mM | 0.3592 mL | 1.7959 mL | 3.5918 mL |
| 10 mM | 0.1796 mL | 0.8979 mL | 1.7959 mL |
| 50 mM | 0.0359 mL | 0.1796 mL | 0.3592 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481