

## 4-hydroxy-3,5-dimethyl-6-[(4E)-4-methylhex-4-en-2-yl]-2-oxooxan-3-yl (2E,4E)-hexa-2,4-dienoate


### Chemical Properties

CAS No. :

Formula: C<sub>20</sub>H<sub>30</sub>O<sub>5</sub>

Molecular Weight: 350.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



### Biological Description

Description	4-hydroxy-3,5-dimethyl-6-[(4E)-4-methylhex-4-en-2-yl]-2-oxooxan-3-yl (2E,4E)-hexa-2,4-dienoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8534 mL	14.267 mL	28.5339 mL
5 mM	0.5707 mL	2.8534 mL	5.7068 mL
10 mM	0.2853 mL	1.4267 mL	2.8534 mL
50 mM	0.0571 mL	0.2853 mL	0.5707 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481