

2-[(19-amino-2,18-dihydroxyicosan-3-yl)oxy]-6-(hydroxymethyl)oxane-3,4,5-triol

Chemical Properties


CAS No. :

Formula: C₂₆H₅₃N₀O₈

Molecular Weight: 507.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[(19-amino-2,18-dihydroxyicosan-3-yl)oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9696 mL	9.8481 mL	19.6963 mL
5 mM	0.3939 mL	1.9696 mL	3.9393 mL
10 mM	0.197 mL	0.9848 mL	1.9696 mL
50 mM	0.0394 mL	0.197 mL	0.3939 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481