

3-hydroxy-2-{7-[3-hydroxy-5-(3-hydroxyheptyl)oxolan-2-yl]heptanamido}propanoic acid

Chemical Properties


CAS No. :

Formula: C₂₁H₃₉N₁O₇

Molecular Weight: 417.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-hydroxy-2-{7-[3-hydroxy-5-(3-hydroxyheptyl)oxolan-2-yl]heptanamido}propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.395 mL	11.9749 mL	23.9498 mL
5 mM	0.479 mL	2.395 mL	4.790 mL
10 mM	0.2395 mL	1.1975 mL	2.395 mL
50 mM	0.0479 mL	0.2395 mL	0.479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481