

2-[(4E,5E)-5-[carboxy(phenyl)methylidene]-1,3-dioxolan-4-ylidene]-2-phenylacetic acid

Chemical Properties


CAS No. :

Formula: C19H14O6

Molecular Weight: 338.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[(4E,5E)-5-[carboxy(phenyl)methylidene]-1,3-dioxolan-4-ylidene]-2-phenylacetic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9558 mL	14.7789 mL	29.5578 mL
5 mM	0.5912 mL	2.9558 mL	5.9116 mL
10 mM	0.2956 mL	1.4779 mL	2.9558 mL
50 mM	0.0591 mL	0.2956 mL	0.5912 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481