

methyl 3,15,23-tris({[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy})tetracosanoate

## Chemical Properties


CAS No. :

Formula: C43H80O20

Molecular Weight: 917.09

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	methyl 3,15,23-tris({[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy})tetracosanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0904 mL	5.452 mL	10.9041 mL
5 mM	0.2181 mL	1.0904 mL	2.1808 mL
10 mM	0.109 mL	0.5452 mL	1.0904 mL
50 mM	0.0218 mL	0.109 mL	0.2181 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481