


3-[(2E,6E)-10,11-dihydroxy-3,7,11-trimethyldodeca-2,6-dien-1-yl]-2,4-dihydroxy-6-methylbenzoic acid

Chemical Properties

CAS No. :

Formula: C₂₃H₃₄O₆

Molecular Weight: 406.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-[(2E,6E)-10,11-dihydroxy-3,7,11-trimethyldodeca-2,6-dien-1-yl]-2,4-dihydroxy-6-methylbenzoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4599 mL	12.2995 mL	24.599 mL
5 mM	0.492 mL	2.4599 mL	4.9198 mL
10 mM	0.246 mL	1.230 mL	2.4599 mL
50 mM	0.0492 mL	0.246 mL	0.492 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481