

m-PEG3-OH

Chemical Properties

CAS No. : 112-35-6

Formula: C7H16O4

Molecular Weight: 164.2

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-PEG3-OH is a PEG-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, or proteolysis-targeting chimeras, are composed of two distinct ligands joined by a linker. One of these ligands binds to a specific protein target, while the other binds to an E3 ubiquitin ligase. When the PROTAC binds to both the target protein and the E3 ligase, it triggers the ubiquitin-proteasome system within cells to degrade the target protein, thereby providing a mechanism for targeted protein degradation.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0901 mL	30.4507 mL	60.9013 mL
5 mM	1.218 mL	6.0901 mL	12.1803 mL
10 mM	0.609 mL	3.0451 mL	6.0901 mL
50 mM	0.1218 mL	0.609 mL	1.218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481