


## 2-Acetyl-4-butyramidophenol

## Chemical Properties

CAS No. :	40188-45-2
Formula:	C <sub>12</sub> H <sub>15</sub> NO <sub>3</sub>
Molecular Weight:	221.25
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Acetyl-4-butyramidophenol is gray crystalline powder.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 55 mg/mL (248.59 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5198 mL	22.5989 mL	45.1977 mL
5 mM	0.904 mL	4.5198 mL	9.0395 mL
10 mM	0.452 mL	2.2599 mL	4.5198 mL
50 mM	0.0904 mL	0.452 mL	0.904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Wiesmann UN, et al. Biochem Biophys Res Commun. 1975 Oct 27;66(4):1338-43.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481