

1-(2,4-Dihydroxyphenyl)-3-(3,4-dihydroxyphenyl)propane, 3'',4'-Dimethylether

Chemical Properties


CAS No. :

Formula: C₁₇H₂₀O₄

Molecular Weight: 288.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

1-(2,4-Dihydroxyphenyl)-3-(3,4-dihydroxyphenyl)propane, 3'',4'-Dimethylether is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4681 mL	17.3406 mL	34.6813 mL
5 mM	0.6936 mL	3.4681 mL	6.9363 mL
10 mM	0.3468 mL	1.7341 mL	3.4681 mL
50 mM	0.0694 mL	0.3468 mL	0.6936 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481