

Cyclocarioside H

Chemical Properties

CAS No. :	1403937-87-0
Formula:	C43H72O13
Molecular Weight:	797.04
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Cyclocarioside H is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2546 mL	6.2732 mL	12.5464 mL
5 mM	0.2509 mL	1.2546 mL	2.5093 mL
10 mM	0.1255 mL	0.6273 mL	1.2546 mL
50 mM	0.0251 mL	0.1255 mL	0.2509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481