


Derrisaponin A

Chemical Properties

CAS No. :	1965272-21-2
Formula:	C56H88O26
Molecular Weight:	1177.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Derrisaponin A is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8494 mL	4.247 mL	8.4941 mL
5 mM	0.1699 mL	0.8494 mL	1.6988 mL
10 mM	0.0849 mL	0.4247 mL	0.8494 mL
50 mM	0.017 mL	0.0849 mL	0.1699 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481