


Deoxytrillenoside A

Chemical Properties

CAS No. :	77658-50-5
Formula:	C47H70O23
Molecular Weight:	1003.05
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Deoxytrillenoside A is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.997 mL	4.9848 mL	9.9696 mL
5 mM	0.1994 mL	0.997 mL	1.9939 mL
10 mM	0.0997 mL	0.4985 mL	0.997 mL
50 mM	0.0199 mL	0.0997 mL	0.1994 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481