


## Lecaniodiscussaponin S4

## Chemical Properties

CAS No. : 84607-60-3  
 Formula: C<sub>46</sub>H<sub>74</sub>O<sub>16</sub>  
 Molecular Weight: 883.08  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Lecaniodiscussaponin S4 is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1324 mL	5.662 mL	11.324 mL
5 mM	0.2265 mL	1.1324 mL	2.2648 mL
10 mM	0.1132 mL	0.5662 mL	1.1324 mL
50 mM	0.0226 mL	0.1132 mL	0.2265 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481