


Abietatriene-3 β -ol

Chemical Properties

CAS No. : 78078-41-8
 Formula: C₂₀H₃₀O
 Molecular Weight: 286.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Abietatriene-3 β -ol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4909 mL	17.4544 mL	34.9089 mL
5 mM	0.6982 mL	3.4909 mL	6.9818 mL
10 mM	0.3491 mL	1.7454 mL	3.4909 mL
50 mM	0.0698 mL	0.3491 mL	0.6982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481