


Bispyrayafoline A

Chemical Properties

CAS No. :	1715114-37-6
Formula:	C ₄₆ H ₄₈ N ₂ O ₄
Molecular Weight:	692.9
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bispyrayafoline A is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4432 mL	7.216 mL	14.4321 mL
5 mM	0.2886 mL	1.4432 mL	2.8864 mL
10 mM	0.1443 mL	0.7216 mL	1.4432 mL
50 mM	0.0289 mL	0.1443 mL	0.2886 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481