


Schirubridilactone D

Chemical Properties

CAS No. : 1210051-95-8
 Formula: C₂₉H₃₄O₁₀
 Molecular Weight: 542.58
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Schirubridilactone D is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.843 mL	9.2152 mL	18.4305 mL
5 mM	0.3686 mL	1.843 mL	3.6861 mL
10 mM	0.1843 mL	0.9215 mL	1.843 mL
50 mM	0.0369 mL	0.1843 mL	0.3686 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481