


## Withametelin B

## Chemical Properties

CAS No. : 123523-09-1  
 Formula: C<sub>28</sub>H<sub>36</sub>O<sub>5</sub>  
 Molecular Weight: 452.59  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Withametelin B is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2095 mL	11.0475 mL	22.0951 mL
5 mM	0.4419 mL	2.2095 mL	4.419 mL
10 mM	0.221 mL	1.1048 mL	2.2095 mL
50 mM	0.0442 mL	0.221 mL	0.4419 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481