


## Plantainoside E

## Chemical Properties

CAS No. :	136083-88-0
Formula:	C <sub>31</sub> H <sub>40</sub> O <sub>16</sub>
Molecular Weight:	668.65
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Plantainoside E is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4956 mL	7.4778 mL	14.9555 mL
5 mM	0.2991 mL	1.4956 mL	2.9911 mL
10 mM	0.1496 mL	0.7478 mL	1.4956 mL
50 mM	0.0299 mL	0.1496 mL	0.2991 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481