


## Isochionoside B

## Chemical Properties

CAS No. :	159116-51-5
Formula:	C36H48O20
Molecular Weight:	800.76
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Isochionoside B is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2488 mL	6.2441 mL	12.4881 mL
5 mM	0.2498 mL	1.2488 mL	2.4976 mL
10 mM	0.1249 mL	0.6244 mL	1.2488 mL
50 mM	0.025 mL	0.1249 mL	0.2498 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481