


## Qinyangshengenin 3-digitoxoside

## Chemical Properties

CAS No. :	1802801-83-7
Formula:	C <sub>34</sub> H <sub>46</sub> O <sub>11</sub>
Molecular Weight:	630.73
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Qinyangshengenin 3-digitoxoside is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5855 mL	7.9273 mL	15.8546 mL
5 mM	0.3171 mL	1.5855 mL	3.1709 mL
10 mM	0.1585 mL	0.7927 mL	1.5855 mL
50 mM	0.0317 mL	0.1585 mL	0.3171 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481