


7-Hydroxykaurenolide

Chemical Properties

CAS No. :	5691-63-4
Formula:	C ₂₀ H ₂₈ O ₃
Molecular Weight:	316.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	7-Hydroxykaurenolide is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1602 mL	15.8008 mL	31.6016 mL
5 mM	0.632 mL	3.1602 mL	6.3203 mL
10 mM	0.316 mL	1.5801 mL	3.1602 mL
50 mM	0.0632 mL	0.316 mL	0.632 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481