

## Compound TN11916

### Chemical Properties


CAS No. :

Formula: C<sub>35</sub>H<sub>46</sub>O<sub>20</sub>

Molecular Weight: 786.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

**Description** (2R,3R,4R,5R,6R)-4-[[[(2S,3R,4R,5R,6S)-4,5-dihydroxy-6-methyl-3-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]oxan-2-yl]oxy]-6-[2-(3,4-dihydroxyphenyl)ethoxy]-5-hydroxy-2-(hydroxymethyl)oxan-3-yl (2E)-3-(3,4-dihydroxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2711 mL	6.3554 mL	12.7108 mL
5 mM	0.2542 mL	1.2711 mL	2.5422 mL
10 mM	0.1271 mL	0.6355 mL	1.2711 mL
50 mM	0.0254 mL	0.1271 mL	0.2542 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481