

2-(3,4-dimethoxyphenoxy)-6-[[[(3,4,5-trihydroxyoxan-2-yl)oxy]methyl]oxane-3,4,5-triol

Chemical Properties


CAS No. :

Formula: C₁₉H₂₈O₁₂

Molecular Weight: 448.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(3,4-dimethoxyphenoxy)-6-[[[(3,4,5-trihydroxyoxan-2-yl)oxy]methyl]oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2301 mL	11.1503 mL	22.3005 mL
5 mM	0.446 mL	2.2301 mL	4.4601 mL
10 mM	0.223 mL	1.115 mL	2.2301 mL
50 mM	0.0446 mL	0.223 mL	0.446 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481