

## Compound TN11927

### Chemical Properties


CAS No. :

Formula: C<sub>24</sub>H<sub>32</sub>O<sub>11</sub>

Molecular Weight: 496.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	(2S,3R,4S,5R)-2-(4-{1,3-dihydroxy-2-[2-hydroxy-4-(3-hydroxypropyl)phenoxy]propyl}-2-methoxyphenoxy)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0141 mL	10.0703 mL	20.1406 mL
5 mM	0.4028 mL	2.0141 mL	4.0281 mL
10 mM	0.2014 mL	1.007 mL	2.0141 mL
50 mM	0.0403 mL	0.2014 mL	0.4028 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481