

2-{4-hydroxy-2,3-bis[(4-hydroxy-3-methoxyphenyl)methyl]butoxy}-6-(hydroxymethyl)oxane-3,4,5-triol


Chemical Properties

CAS No. :

Formula: C₂₆H₃₆O₁₁

Molecular Weight: 524.56

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-{4-hydroxy-2,3-bis[(4-hydroxy-3-methoxyphenyl)methyl]butoxy}-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9064 mL	9.5318 mL	19.0636 mL
5 mM	0.3813 mL	1.9064 mL	3.8127 mL
10 mM	0.1906 mL	0.9532 mL	1.9064 mL
50 mM	0.0381 mL	0.1906 mL	0.3813 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481