

Compound TN11956

Chemical Properties


CAS No. :

Formula: C₂₃H₃₆O₅

Molecular Weight: 392.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(2E)-6-(5-hydroxy-4-((2-methylpropanoyl)oxy)methyl)cyclohex-3-en-1-yl)-2-methylhepta-2,6-dien-1-yl 2-methylpropanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5475 mL	12.7376 mL	25.4751 mL
5 mM	0.5095 mL	2.5475 mL	5.095 mL
10 mM	0.2548 mL	1.2738 mL	2.5475 mL
50 mM	0.051 mL	0.2548 mL	0.5095 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481