

## 5,7-dihydroxy-2-phenyl-8-[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]-4H-chromen-4-one

### Chemical Properties


CAS No. :

Formula: C<sub>21</sub>H<sub>20</sub>O<sub>9</sub>

Molecular Weight: 416.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	5,7-dihydroxy-2-phenyl-8-[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]-4H-chromen-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4017 mL	12.0083 mL	24.0165 mL
5 mM	0.4803 mL	2.4017 mL	4.8033 mL
10 mM	0.2402 mL	1.2008 mL	2.4017 mL
50 mM	0.048 mL	0.2402 mL	0.4803 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481