

(1S)-1,7,7-trimethylbicyclo[2.2.1]heptan-2-yl (2E)-3-(4-hydroxyphenyl)prop-2-enoate

Chemical Properties


CAS No. :

Formula: C₁₉H₂₄O₃

Molecular Weight: 300.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(1S)-1,7,7-trimethylbicyclo[2.2.1]heptan-2-yl (2E)-3-(4-hydroxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3289 mL	16.6445 mL	33.2889 mL
5 mM	0.6658 mL	3.3289 mL	6.6578 mL
10 mM	0.3329 mL	1.6644 mL	3.3289 mL
50 mM	0.0666 mL	0.3329 mL	0.6658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481