

Compound TN11993

Chemical Properties


CAS No. :

Formula: C63H102O31

Molecular Weight: 1355.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-[[[(2S,3R,4S,5R,6S)-5-[[[(2R,3R,4R)-3,4-dihydroxy-4-(hydroxymethyl)oxolan-2-yl]oxy]-3-hydroxy-4-[(3,4,5-trihydroxy-6-methyloxan-2-yl)oxy]oxan-2-yl]oxy]-3,4-dihydroxy-6-methyloxan-2-yl]oxy]-4,5-dihydroxyoxan-2-yl (4aS,6aS,6bR,8R,9R,10R,11S,12aR)-8,11-dihydroxy-9-(hydroxymethyl)-2,2,6a,6b,9,12a-hexamethyl-10-[[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7377 mL	3.6887 mL	7.3775 mL
5 mM	0.1475 mL	0.7377 mL	1.4755 mL
10 mM	0.0738 mL	0.3689 mL	0.7377 mL
50 mM	0.0148 mL	0.0738 mL	0.1475 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481