

Compound TN12014

Chemical Properties


CAS No. :

Formula: C₂₂H₃₆O₈

Molecular Weight: 428.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-[2-(3,4-dihydroxy-2,5-dimethoxyoxolan-3-yl)ethyl]-8-hydroxy-5,6,8a-trimethyl-3,4,4a,5,6,7,8,8a-octahydronaphthalene-1-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3336 mL	11.6681 mL	23.3361 mL
5 mM	0.4667 mL	2.3336 mL	4.6672 mL
10 mM	0.2334 mL	1.1668 mL	2.3336 mL
50 mM	0.0467 mL	0.2334 mL	0.4667 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481