

Compound TN12021

Chemical Properties


CAS No. :

Formula: C54H86O23

Molecular Weight: 1103.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3S,4S,5R,6R)-6-(((3S,6aR,6bS,8aS,14bR)-4,4,6a,6b,11,11,14b-heptamethyl-8a-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)carbonyl)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy)-4-(((2S,3R,4R,5R,6S)-4,5-dihydroxy-6-methyl-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)oxan-2-yl]oxy)-3,5-dihydroxyoxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9064 mL	4.532 mL	9.064 mL
5 mM	0.1813 mL	0.9064 mL	1.8128 mL
10 mM	0.0906 mL	0.4532 mL	0.9064 mL
50 mM	0.0181 mL	0.0906 mL	0.1813 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481